

Mouse Adar/Double-stranded RNA-specific adenosine deaminase ELISA Kit



Catalog No: #EK19451

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Adar/Double-stranded RNA-specific adenosine deaminase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Double-stranded RNA-specific adenosine deaminase
Target Name	Double-stranded RNA-specific adenosine deaminase
Other Names	DRADA, RNA adenosine deaminase 1
Accession No.	Q99MU3
Uniprot	Q99MU3
GeneID	56417;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.069ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Adar

Note: This product is for in vitro research use only