

Rat Bgn/Biglycan ELISA Kit

Catalog No: #EK19452

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Bgn/Biglycan ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Biglycan
Target Name	Biglycan
Other Names	Bone/cartilage proteoglycan I,PG-S1
Accession No.	P47853
Uniprot	P47853
GeneID	25181;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.086ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Bgn

Note: This product is for in vitro research use only