

Human CELA2A/Chymotrypsin-like elastase family member 2A ELISA Kit



Catalog No: #EK19453

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Human CELA2A/Chymotrypsin-like elastase family member 2A ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Hu |
| Specificity | Natural and recombinant Human Chymotrypsin-like elastase family member 2A |
| Target Name | Chymotrypsin-like elastase family member 2A |
| Other Names | Elastase-2A,ELA2A |
| Accession No. | P08217 |
| Uniprot | P08217 |
| GeneID | 63036; |

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.403ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CELA2A

Note: This product is for in vitro research use only