

# Mouse Birc2/Baculoviral IAP repeat-containing protein 2 ELISA Kit



Catalog No: #EK19455

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Birc2/Baculoviral IAP repeat-containing protein 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Baculoviral IAP repeat-containing protein 2
Target Name	Baculoviral IAP repeat-containing protein 2
Other Names	Cellular inhibitor of apoptosis 1,C-IAP1,Inhibitor of apoptosis protein 2,miAP2,RING-type E3 ubiquitin transferase BIRC2
Accession No.	Q62210
Uniprot	Q62210
GeneID	11797;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.12ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Birc2

Note: This product is for in vitro research use only