

# Mouse Bco2/Beta,beta-carotene 9',10'-oxygenase ELISA Kit



Catalog No: #EK19456

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Bco2/Beta,beta-carotene 9',10'-oxygenase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Beta,beta-carotene 9',10'-oxygenase
Target Name	Beta,beta-carotene 9',10'-oxygenase
Other Names	B-diox-II,Beta-carotene dioxygenase 2,Bcdo2
Accession No.	Q99NF1
Uniprot	Q99NF1
GeneID	170752;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.11ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Bco2

Note: This product is for in vitro research use only