

# Bovine ALPL/Alkaline phosphatase, tissue-nonspecific isozyme ELISA Kit



Catalog No: #EK19457

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine ALPL/Alkaline phosphatase, tissue-nonspecific isozyme ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Alkaline phosphatase, tissue-nonspecific isozyme
Target Name	Alkaline phosphatase, tissue-nonspecific isozyme
Other Names	AP-TNAP,Alkaline phosphatase liver/bone/kidney isozyme
Accession No.	P09487
Uniprot	P09487
GeneID	280994;

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:ALPL

Note: This product is for in vitro research use only