

# Human ADAD1/Adenosine deaminase domain-containing protein 1 ELISA Kit

Catalog No: #EK19459

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ADAD1/Adenosine deaminase domain-containing protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Adenosine deaminase domain-containing protein 1
Target Name	Adenosine deaminase domain-containing protein 1
Other Names	Testis nuclear RNA-binding protein, TENR
Accession No.	Q96M93
Uniprot	Q96M93
GeneID	132612;

## Application Details

Assay Type: Sandwich Detect Range: 6.25-400pg/mL Sensitivity: 2.1pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: ADAD1

Note: This product is for in vitro research use only