

Mouse Adad1/Adenosine deaminase domain-containing protein 1 ELISA Kit

Catalog No: #EK19460

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Adad1/Adenosine deaminase domain-containing protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Adenosine deaminase domain-containing protein 1
Target Name	Adenosine deaminase domain-containing protein 1
Other Names	Testis nuclear RNA-binding protein, Tenr
Accession No.	Q5SUE7
Uniprot	Q5SUE7
GeneID	21744;

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 41pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Adad1

Note: This product is for in vitro research use only