

Rat Adad1/Adenosine deaminase domain-containing protein 1 ELISA Kit



Catalog No: #EK19461

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Adad1/Adenosine deaminase domain-containing protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Adenosine deaminase domain-containing protein 1
Target Name	Adenosine deaminase domain-containing protein 1
Other Names	Testis nuclear RNA-binding protein, Tenr
Accession No.	Q3KR54
Uniprot	Q3KR54
GeneID	294979;

Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.789ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Adad1

Note: This product is for in vitro research use only