

Human CALB1/Calbindin ELISA Kit

Catalog No: #EK19462

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CALB1/Calbindin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Calbindin
Target Name	Calbindin
Other Names	Calbindin D28,D-28K,Vitamin D-dependent calcium-binding protein, avian-type,CAB27
Accession No.	P05937
Uniprot	P05937
GeneID	793;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.111ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CALB1

Note: This product is for in vitro research use only