

Mouse Cadm3/Cell adhesion molecule 3 ELISA Kit

Catalog No: #EK19464

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Mouse Cadm3/Cell adhesion molecule 3 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Ms |
| Specificity | Natural and recombinant Mouse Cell adhesion molecule 3 |
| Target Name | Cell adhesion molecule 3 |
| Other Names | Immunoglobulin superfamily member 4B,IgSF4B,Nectin-like protein 1,NECL-1,Synaptic cell adhesion molecule 3,TSLC1-like protein 1,Igsf4b,Necl1,Syncam3,Tsll1 |
| Accession No. | Q99N28 |
| Uniprot | Q99N28 |
| GeneID | 94332; |

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.21ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cadm3

Note: This product is for in vitro research use only