

Rat Cadm3/Cell adhesion molecule 3 ELISA Kit

Catalog No: #EK19465

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cadm3/Cell adhesion molecule 3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cell adhesion molecule 3
Target Name	Cell adhesion molecule 3
Other Names	Immunoglobulin superfamily member 4B,IgSF4B,Nectin-like protein 1,NECL-1,Igsf4b,Necl1
Accession No.	Q1WIM3
Uniprot	Q1WIM3
GeneID	360882;

Application Details

Assay Type: Sandwich Detect Range:1.56-100ng/mL Sensitivity:0.86ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cadm3

Note: This product is for in vitro research use only