

## Mouse Ca4/Carbonic anhydrase 4 ELISA Kit

Catalog No: #EK19466

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Ca4/Carbonic anhydrase 4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Carbonic anhydrase 4
Target Name	Carbonic anhydrase 4
Other Names	Carbonate dehydratase IV,Carbonic anhydrase IV,CA-IV,Car4
Accession No.	Q64444
Uniprot	Q64444
GeneID	12351;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.081ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Ca4

Note: This product is for in vitro research use only