

Rat Ca4/Carbonic anhydrase 4 ELISA Kit

Catalog No: #EK19467

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Ca4/Carbonic anhydrase 4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Carbonic anhydrase 4
Target Name	Carbonic anhydrase 4
Other Names	Carbonate dehydratase IV,Carbonic anhydrase IV,CA-IV
Accession No.	P48284
Uniprot	P48284
GeneID	29242;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.12ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Ca4

Note: This product is for in vitro research use only