

# Mouse Bcl6b/B-cell CLL/lymphoma 6 member B protein ELISA Kit

Catalog No: #EK19468

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Bcl6b/B-cell CLL/lymphoma 6 member B protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse B-cell CLL/lymphoma 6 member B protein
Target Name	B-cell CLL/lymphoma 6 member B protein
Other Names	Bcl6-associated zinc finger protein,Bazf
Accession No.	O88282
Uniprot	O88282
GeneID	12029;

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.084ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Bcl6b

Note: This product is for in vitro research use only