

Mouse Alk/ALK tyrosine kinase receptor ELISA Kit

Catalog No: #EK19471

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Alk/ALK tyrosine kinase receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse ALK tyrosine kinase receptor
Target Name	ALK tyrosine kinase receptor
Other Names	Anaplastic lymphoma kinase,CD246
Accession No.	P97793
Uniprot	P97793
GeneID	11682;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:41.4pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Alk

Note: This product is for in vitro research use only