

Mouse Cdk16/Cyclin-dependent kinase 16 ELISA Kit

Catalog No: #EK19473



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Mouse Cdk16/Cyclin-dependent kinase 16 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cyclin-dependent kinase 16
Target Name	Cyclin-dependent kinase 16
Other Names	CRK5,Cell division protein kinase 16,PCTAIRE-motif protein kinase 1,Serine/threonine-protein kinase PCTAIRE-1,Crk5,Pctaire1,Pctk1
Accession No.	Q04735
Uniprot	Q04735
GeneID	18555;

Application Details

Assay Type: SandwichDetect Range:1.56-100ng/mL Sensitivity:0.46ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Cdk16

Note: This product is for in vitro research use only