

Mouse Bmp10/Bone morphogenetic protein 10 ELISA Kit



Catalog No: #EK19474

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Bmp10/Bone morphogenetic protein 10 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Bone morphogenetic protein 10
Target Name	Bone morphogenetic protein 10
Other Names	BMP-10
Accession No.	Q9R229
Uniprot	Q9R229
GeneID	12154;

Application Details

Assay Type: Sandwich Detect Range: 2.5-160ng/mL Sensitivity: 1.25ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Bmp10

Note: This product is for in vitro research use only