

Human CENPB/Major centromere autoantigen B ELISA Kit



Catalog No: #EK19478

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CENPB/Major centromere autoantigen B ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Major centromere autoantigen B
Target Name	Major centromere autoantigen B
Other Names	Centromere protein B,CENP-B
Accession No.	P07199
Uniprot	P07199
GeneID	1059;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.12ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CENPB

Note: This product is for in vitro research use only