

Mouse Atp1a3/Sodium/potassium-transporting ATPase subunit alpha-3 ELISA Kit



Catalog No: #EK19479

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Atp1a3/Sodium/potassium-transporting ATPase subunit alpha-3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Sodium/potassium-transporting ATPase subunit alpha-3
Target Name	Sodium/potassium-transporting ATPase subunit alpha-3
Other Names	Na(+)/K(+) ATPase alpha-3 subunit,Na(+)/K(+) ATPase alpha(III) subunit,Sodium pump subunit alpha-3
Accession No.	Q6PIC6
Uniprot	Q6PIC6
GeneID	232975;

Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Atp1a3

Note: This product is for in vitro research use only