

Human ANXA10/Annexin A10 ELISA Kit

Catalog No: #EK19480

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ANXA10/Annexin A10 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Annexin A10
Target Name	Annexin A10
Other Names	Annexin-10,Annexin-14,ANX14
Accession No.	Q9UJ72
Uniprot	Q9UJ72
GeneID	11199;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:ANXA10

Note: This product is for in vitro research use only