

## Human CD52/CAMPATH-1 antigen ELISA Kit

Catalog No: #EK19481

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CD52/CAMPATH-1 antigen ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human CAMPATH-1 antigen
Target Name	CAMPATH-1 antigen
Other Names	CDw52,Cambridge pathology 1 antigen,Epididymal secretory protein E5,Human epididymis-specific protein 5,He5,CD52,CDW52,HE5
Accession No.	P31358
Uniprot	P31358
GeneID	1043;

## Application Details

Assay Type: SandwichDetect Range:0.156-10ng/ml Sensitivity:0.082ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:CD52

Note: This product is for in vitro research use only