

Mouse Cd52/CAMPATH-1 antigen ELISA Kit

Catalog No: #EK19482



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Mouse Cd52/CAMPATH-1 antigen ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse CAMPATH-1 antigen
Target Name	CAMPATH-1 antigen
Other Names	Lymphocyte differentiation antigen B7,CD52,Cdw52,Mb7
Accession No.	Q64389
Uniprot	Q64389
GeneID	23833;

Application Details

Assay Type: SandwichDetect Range:0.156-10ng/ml Sensitivity:0.083ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Cd52

Note: This product is for in vitro research use only