

## Rat Cd52/CAMPATH-1 antigen ELISA Kit

Catalog No: #EK19483

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Cd52/CAMPATH-1 antigen ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat CAMPATH-1 antigen
Target Name	CAMPATH-1 antigen
Other Names	Lymphocyte differentiation antigen B7, CD52, Cdw52, Rb7
Accession No.	Q63064
Uniprot	Q63064
GeneID	117054;

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.112ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cd52

Note: This product is for in vitro research use only