

Mouse Aldh5a1/Succinate-semialdehyde dehydrogenase, mitochondrial ELISA Kit

Catalog No: #EK19484

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Aldh5a1/Succinate-semialdehyde dehydrogenase, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Succinate-semialdehyde dehydrogenase, mitochondrial
Target Name	Succinate-semialdehyde dehydrogenase, mitochondrial
Other Names	Aldehyde dehydrogenase family 5 member A1,NAD(+)-dependent succinic semialdehyde dehydrogenase
Accession No.	Q8BWF0
Uniprot	Q8BWF0
GeneID	214579;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.154ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Aldh5a1

Note: This product is for in vitro research use only