

## Human ADH1C/Alcohol dehydrogenase 1C ELISA Kit

Catalog No: #EK19487

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ADH1C/Alcohol dehydrogenase 1C ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Alcohol dehydrogenase 1C
Target Name	Alcohol dehydrogenase 1C
Other Names	Alcohol dehydrogenase subunit gamma,ADH3
Accession No.	P00326
Uniprot	P00326
GeneID	126;

## Application Details

Assay Type: Sandwich Detect Range:15.6-1000pg/mL Sensitivity:7.4pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:ADH1C

Note: This product is for in vitro research use only