

# Mouse Chk1/Serine/threonine-protein kinase Chk1 ELISA Kit



Catalog No: #EK19488

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Chk1/Serine/threonine-protein kinase Chk1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Serine/threonine-protein kinase Chk1
Target Name	Serine/threonine-protein kinase Chk1
Other Names	CHK1 checkpoint homolog,Checkpoint kinase-1,Chk1
Accession No.	O35280
Uniprot	O35280
GeneID	12649;

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 43pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Chk1

Note: This product is for in vitro research use only