

## Mouse Ago2/Protein argonaute-2 ELISA Kit

Catalog No: #EK19489

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Ago2/Protein argonaute-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Protein argonaute-2
Target Name	Protein argonaute-2
Other Names	Argonaute2, Argonaute RISC catalytic component 2, Eukaryotic translation initiation factor 2C 2, eIF-2C 2, Piwi/argonaute family protein melf2c2, Protein slicer, Eif2c2, Kiaa4215
Accession No.	Q8CJG0
Uniprot	Q8CJG0
GeneID	239528;

## Application Details

Assay Type: Sandwich Detect Range: 78-5000 pg/mL Sensitivity: 39.4 pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Ago2

Note: This product is for in vitro research use only