

## Rat Ago2/Protein argonaute-2 ELISA Kit

Catalog No: #EK19490

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Ago2/Protein argonaute-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Protein argonaute-2
Target Name	Protein argonaute-2
Other Names	Argonaute2, Argonaute RISC catalytic component 2, Eukaryotic translation initiation factor 2C 2, eIF-2C 2, Golgi ER protein 95 kDa, GERp95, Protein slicer, Eif2c2
Accession No.	Q9QZ81
Uniprot	Q9QZ81
GeneID	59117;

## Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.174ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Ago2

Note: This product is for in vitro research use only