

## Human CD69/Early activation antigen CD69 ELISA Kit

Catalog No: #EK19491

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CD69/Early activation antigen CD69 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Early activation antigen CD69
Target Name	Early activation antigen CD69
Other Names	Activation inducer molecule,AIM,BL-AC/P26,C-type lectin domain family 2 member C,EA1,Early T-cell activation antigen p60,GP32/28,Leukocyte surface antigen Leu-23,MLR-3,CD69,CLEC2C
Accession No.	Q07108
Uniprot	Q07108
GeneID	969;

## Application Details

Assay Type: SandwichDetect Range:0.312-20ng/mL Sensitivity:0.153ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene\_full\_name:CD69

Note: This product is for in vitro research use only