

Human ANTXR2/Anthrax toxin receptor 2 ELISA Kit

Catalog No: #EK19494

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human ANTXR2/Anthrax toxin receptor 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Anthrax toxin receptor 2
Target Name	Anthrax toxin receptor 2
Other Names	Capillary morphogenesis gene 2 protein, CMG-2, CMG2
Accession No.	P58335
Uniprot	P58335
GeneID	118429;

Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/mL Sensitivity: 0.78ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: ANTXR2

Note: This product is for in vitro research use only