

Mouse Bai3/Brain-specific angiogenesis inhibitor 3 ELISA Kit

Catalog No: #EK19495

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Bai3/Brain-specific angiogenesis inhibitor 3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Brain-specific angiogenesis inhibitor 3
Target Name	Brain-specific angiogenesis inhibitor 3
Accession No.	Q80ZF8
Uniprot	Q80ZF8
GeneID	210933;

Application Details

Assay Type: SandwichDetect Range:78-5000pg/mL Sensitivity:39 pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Bai3

Note: This product is for in vitro research use only