

Mouse Cldn2/Claudin-2 ELISA Kit

Catalog No: #EK19497



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Cldn2/Claudin-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Claudin-2
Target Name	Claudin-2
Accession No.	O88552
Uniprot	O88552
GeneID	12738;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.084ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cldn2

Note: This product is for in vitro research use only