

# Bovine CLDN3/Claudin-3 ELISA Kit

Catalog No: #EK19498



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | Bovine CLDN3/Claudin-3 ELISA Kit         |
| Brief Description  | ELISA Kits                               |
| Applications       | ELISA                                    |
| Species Reactivity | Bovine                                   |
| Specificity        | Natural and recombinant Bovine Claudin-3 |
| Target Name        | Claudin-3                                |
| Accession No.      | Q765N9                                   |
| Uniprot            | Q765N9                                   |
| GeneID             | 404153;                                  |

## Application Details

Assay Type: Sandwich Detect Range: 0.625-40ng/mL Sensitivity: 0.38ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CLDN3

Note: This product is for in vitro research use only