

Mouse Cldn3/Claudin-3 ELISA Kit

Catalog No: #EK19499

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cldn3/Claudin-3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Claudin-3
Target Name	Claudin-3
Other Names	Clostridium perfringens enterotoxin receptor 2,CPE-R 2,Cpetr2
Accession No.	Q9Z0G9
Uniprot	Q9Z0G9
GeneID	12739;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.085ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cldn3

Note: This product is for in vitro research use only