

TAS2R44 Antibody

Catalog No: #44658

Package Size: #44658-1 50ul #44658-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TAS2R44 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	TAS2R44 Antibody detects endogenous levels of total TAS2R44
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human TAS2R44
Target Name	TAS2R44
Other Names	TAS2R31, T2R53, TAS2R44, T2R31, T2R44
Accession No.	Swiss-Prot#: P59538NCBI Gene ID:259290
Uniprot	P59538
GeneID	259290;
Calculated MW	35 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only