

Human CPB1/Carboxypeptidase B ELISA Kit

Catalog No: #EK19503

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CPB1/Carboxypeptidase B ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Carboxypeptidase B
Target Name	Carboxypeptidase B
Other Names	Pancreas-specific protein,PASP,CPB,PCPB
Accession No.	P15086
Uniprot	P15086
GeneID	1360;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.069ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CPB1

Note: This product is for in vitro research use only