

## Mouse Crygd/Gamma-crystallin D ELISA Kit

Catalog No: #EK19504

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Crygd/Gamma-crystallin D ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Gamma-crystallin D
Target Name	Gamma-crystallin D
Other Names	Gamma-D-crystallin,Gamma-crystallin 1
Accession No.	P04342
Uniprot	P04342
GeneID	12967;

## Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.36ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Crygd

Note: This product is for in vitro research use only