

Mouse Crym/Ketimine reductase mu-crystallin ELISA Kit



Catalog No: #EK19505

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Crym/Ketimine reductase mu-crystallin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Ketimine reductase mu-crystallin
Target Name	Ketimine reductase mu-crystallin
Other Names	NADP-regulated thyroid-hormone-binding protein
Accession No.	O54983
Uniprot	O54983
GeneID	12971;

Application Details

Assay Type: Sandwich Detect Range: 1.56-100ng/ml Sensitivity: 0.782ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Crym

Note: This product is for in vitro research use only