

Rat Crym/Ketimine reductase mu-crystallin ELISA Kit

Catalog No: #EK19506

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Crym/Ketimine reductase mu-crystallin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Ketimine reductase mu-crystallin
Target Name	Ketimine reductase mu-crystallin
Other Names	CDK108,NADP-regulated thyroid-hormone-binding protein
Accession No.	Q9QYU4
Uniprot	Q9QYU4
GeneID	117024;

Application Details

Assay Type: Sandwich Detect Range:15.6-1000pg/mL Sensitivity:5.4pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Crym

Note: This product is for in vitro research use only