

Human CALCRL/Calcitonin gene-related peptide type 1 receptor ELISA Kit



Catalog No: #EK19507

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CALCRL/Calcitonin gene-related peptide type 1 receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Calcitonin gene-related peptide type 1 receptor
Target Name	Calcitonin gene-related peptide type 1 receptor
Other Names	CGRP type 1 receptor, Calcitonin receptor-like receptor, CGRPR
Accession No.	Q16602
Uniprot	Q16602
GeneID	10203;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.057ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CALCRL

Note: This product is for in vitro research use only