

Mouse Clu/Clusterin ELISA Kit

Catalog No: #EK19508

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Clu/Clusterin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Clusterin
Target Name	Clusterin
Other Names	Apolipoprotein J,Apo-J,Clustrin,Sulfated glycoprotein 2,SGP-2,Apoj,Msgp-2
Accession No.	Q06890
Uniprot	Q06890
GeneID	12759;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.063ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Clu

Note: This product is for in vitro research use only