

## Pig CLU/Clusterin ELISA Kit

Catalog No: #EK19509

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Pig CLU/Clusterin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Pig
Specificity	Natural and recombinant Pig Clusterin
Target Name	Clusterin
Other Names	CP40, Complement cytolysis inhibitor, CLI
Accession No.	Q29549
Uniprot	Q29549
GeneID	397025;

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.059ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CLU

Note: This product is for in vitro research use only