

Mouse Abcb9/ATP-binding cassette sub-family B member 9 ELISA Kit



Catalog No: #EK19510

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Abcb9/ATP-binding cassette sub-family B member 9 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse ATP-binding cassette sub-family B member 9
Target Name	ATP-binding cassette sub-family B member 9
Other Names	ATP-binding cassette transporter 9,ABC transporter 9 protein,TAP-like protein,TAPL,Kiaa1520
Accession No.	Q9JJ59
Uniprot	Q9JJ59
GeneID	56325;

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.429ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Abcb9

Note: This product is for in vitro research use only