

Mouse Cndp2/Cytosolic non-specific dipeptidase ELISA Kit



Catalog No: #EK19513

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cndp2/Cytosolic non-specific dipeptidase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cytosolic non-specific dipeptidase
Target Name	Cytosolic non-specific dipeptidase
Other Names	CNDP dipeptidase 2, Glutamate carboxypeptidase-like protein 1, Cn2
Accession No.	Q9D1A2
Uniprot	Q9D1A2
GeneID	66054;

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.401ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cndp2

Note: This product is for in vitro research use only