

Rat Cndp2/Cytosolic non-specific dipeptidase ELISA Kit



Catalog No: #EK19514

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cndp2/Cytosolic non-specific dipeptidase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Cytosolic non-specific dipeptidase
Target Name	Cytosolic non-specific dipeptidase
Other Names	CNDP dipeptidase 2
Accession No.	Q6Q0N1
Uniprot	Q6Q0N1
GeneID	291394;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.061ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cndp2

Note: This product is for in vitro research use only