

Mouse Cndp1/Beta-Ala-His dipeptidase ELISA Kit

Catalog No: #EK19515

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cndp1/Beta-Ala-His dipeptidase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Beta-Ala-His dipeptidase
Target Name	Beta-Ala-His dipeptidase
Other Names	CNDP dipeptidase 1,Carnosine dipeptidase 1,Cn1
Accession No.	Q8BUG2
Uniprot	Q8BUG2
GeneID	338403;

Application Details

Assay Type: Sandwich Detect Range:0.78-50ng/mL Sensitivity:0.34ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cndp1

Note: This product is for in vitro research use only