

Rat Cndp1/Beta-Ala-His dipeptidase ELISA Kit

Catalog No: #EK19516

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Cndp1/Beta-Ala-His dipeptidase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Beta-Ala-His dipeptidase
Target Name	Beta-Ala-His dipeptidase
Other Names	CNDP dipeptidase 1, Carnosine dipeptidase 1, Cn1
Accession No.	Q66HG3
Uniprot	Q66HG3
GeneID	307212;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.081ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cndp1

Note: This product is for in vitro research use only