

# Bovine CSH1/Chorionic somatomammotropin hormone 1 ELISA Kit



Catalog No: #EK19517

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine CSH1/Chorionic somatomammotropin hormone 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Chorionic somatomammotropin hormone 1
Target Name	Chorionic somatomammotropin hormone 1
Other Names	BPLP-I,Placental lactogen I,PL1
Accession No.	P09611
Uniprot	P09611

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.157ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CSH1

Note: This product is for in vitro research use only