

Bovine ACHE/Acetylcholinesterase ELISA Kit

Catalog No: #EK19518

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine ACHE/Acetylcholinesterase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Acetylcholinesterase
Target Name	Acetylcholinesterase
Other Names	AChE
Accession No.	P23795
Uniprot	P23795
GeneID	540446;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.158ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:ACHE

Note: This product is for in vitro research use only